


Externally Powered Systems



The NY Electric Elbow fills the need for a rugged, reliable, and moderately priced externally-powered elbow that can be used by children as well as adults. It has been clinically tested and proven in applications around the world and offers significant benefits, particularly to the bilateral Shoulder Disarticulation or Above Elbow patient.

The electric elbow, available in exo or endoskeletal, is compatible with the Synergetic Prehensor, the Prehension Actuator, and the Michigan Electric Hook. An interchangeable variety of switching systems and myoelectric controls can be used individually or in combination to operate the externally powered components. Commonly available, rechargeable, environmentally safe, Nickel Metal Hydride batteries will power a device for a full day under normal operating conditions.

Fabrication of a complete prosthesis is also available through the Central Fabrication Department of Hosmer Dorrance.




Externally Powered Systems

5V Battery Pack w/Batteries.....	H115
5V Remote Battery Pack.....	H115
6V Modular Battery Pack.....	H115
Battery Chargers	H116
Control Systems.....	H106
Dual-in-Line 5V Remote Battery Pack	H115
Electric Accessories	H117
Fabrication	H112
Li-Ion Power System Battery	H116
Michigan Electric Hook.....	H110
Myoelectric Accessories	H117
Myoelectric Control Package	H109
Myoelectric Control Processor Assembly	H112
MyoLab II Analyzer	H117
NU-VA Synergetic Prehensor Specifications	H113
NU-VA Synergetic Prehensor Terminal Device	H111
NY Electric Elbow.....	H105
NY Electric Elbow Specifications	H107
NY Electric Elbow Switch Controls.....	H114
NY Prehension Actuator (P.A.)	H108
Prehension Actuator or Michigan Hook Myoelectric Control.....	H114
Prehension Actuator or Michigan Hook Switch Controls.....	H114
Replacement Batteries.....	H116
Switch Controls.....	H109
Wrists with Electrical Connectors	H112

NY Electric Elbow

NY Electric Elbow



- For upper extremity amputees with special requirements, such as limited excursion
- May be controlled myoelectrically or with a switch

Features:

- Modular designs of elbow and accessories allow many prostheses configurations
- Interchangeable with the E-200 and E-400 style conventional elbows
- Exo or Endoskeletal application
- Compatible with the NY Prehension Actuator, the Michigan Child Electric Hook and the NU-VA Synergetic Prehensor
- Myoelectric or switch controlled
- Free swing at full extension
- Lightweight; extremely quiet
- Low maintenance
- Commonly available rechargeable NiMH batteries operate the elbow for a full day

NY Electric Elbow, Exoskeletal

- 57431 Exoskeletal Electric Elbow, Medium, Free Swing
- 57432 Exoskeletal Electric Elbow, Large, Free Swing



NY Electric Elbow, Endoskeletal

- 57442 Endoskeletal Electric Elbow, Medium, Free Swing
- 57443 Endoskeletal Electric Elbow, Large, Free Swing



NY Electric Elbow System Components

Control Systems



Switch Controls

For NY Electric Elbows

- 57412 Two Position Pull Switch, NY
- 71987 Two Position/Two Button Push Switch, OCCC
- 56277 Two Position/Two Button Push Switch, NY
- 71977 Three Function Pull Switch, NY (Not Shown)
- 58243 Switch, On/Off Push Pad



Modular Battery Pack

Batteries included

- 56891 5V, Remote Battery Pack
- 56356 5V, Remote Battery Pack, Dual-In-Line
- 54729 5V, Battery Pack in Pouch
- 56776 6V, Endoskeletal Battery Pack
- 61859 7.2V, Li-Ion Battery Pack

Includes:

- A) 61589 Charge Adapter
- B) 61584 7.2v, Li-Ion Battery

Battery Chargers

- 54726 5V Battery Charger, 115V
- 56993 5V or 6V Battery Charger, 115V
- 71617 5V or 6V Battery Charger, 220V (For International Use)
- 60949 Timer for Battery Charger
- 61858 7.2V, Li-Ion Battery Charger Set


Optional

Please Call Central Fabrication For Information


- 58200 Myo Analyzer Kit
- 57818 Laminating Form for EMG Electrode
- 57820 Laminating Form for Ground Electrode
- 57816 Electrode Button, Low Profile, .09" High
- 57815 Electrode Button, Medium Profile, .19" High
- 57814 Electrode Button, High Profile, .28" High
- 57433 Retrofit Electric Elbow w/Free Swing

NY Electric Elbow Specifications

NY Electric Elbow Specifications



Elbow Torque Under Forearm Load	Lift Time (flexion) 5° - 135°	Current	Free Forearm Swing
0 in-lbs	1.3 seconds	.36 amps	
15 in-lbs	2.3 seconds	1.1 amps	
28 in-lbs	4.9 seconds	1.7 amps	5° - 135°
30 in-lbs	STALL	2.0 amps	
Break Away Clutch 18-20 ft-lbs	OVERLOAD	0	



Power Requirements:

Type: AA NiMH battery cell, 1.3 volts each @ 1200mAh, Alkaline batteries can be used for temporary operation but should never be recharged.

Voltage: 5.2 volts = 4 batteries

6.5 volts = 5 batteries

7.2 volts = New 7.2 Li-lone battery pack

Charging Rate: Full charge @ rate of 120mA for 10 hours; charging the batteries can be easily done using the charging jack without removing batteries from prosthesis.


Charger: USE ONLY HOSMER RECOMMENDED CHARGERS WITH THESE POWER PACKS. CHARGING CURRENT, VOLTAGE, AND POLARITY MUST MATCH.

NY Electric Elbow Specifications:

- **Myoelectric Control:** Two EMG sites required, flexion and extension, with ground
- **Switch Control:** Two function control switch, flexion and extension, push or pull switches
- **Weight:** Medium elbow (E-200 size) 15.5oz (439g); large elbow (E-400 size) 16oz, (453g)
- **Height:** Large or medium elbow axis to socket end: 2-1/2" (6.4cm)
- **Diameter:** Socket interface (turntable) (measurement A): medium elbow, 2-3/8" (6cm); large elbow, 2-13/16" (7.1cm)
- **Saddle Width, (measurement B):** Medium elbow, 2-1/4" (5.7cm); large elbow, 2-3/8" (6cm)

NY Prehension Actuator (P.A.)

NY Prehension Actuator (P.A.)




- Motor assisted forearm allows patient with limited excursion to more easily operate a variety of terminal devices
- Significantly decreases the amount of excursion required to open the terminal device
- Forearm custom fabricated to your specifications and requirements

Features:

- Switch or myoelectrically controlled
- Manual pronation/supination with adjustable friction allowed from mid-forearm turntable
- Operates for several days on one charge of commonly available AA NiMH batteries
- Power cutoff feature prevents inadvertent draining of power pack
- Low maintenance
- Available in two wrist sizes (1-3/4" and 2") and compatible with all adult elbows

	Dimension A
NY Elbow, Medium	2.25"(5.72cm)
NY Elbow, Large	2.38"(6.05cm)
E-200 Elbow, Medium	2.12"(5.38cm)
E-400 Elbow, Large	2.38"(6.05cm)



P.A. Configurations for Fabrication

- Wrist and elbow saddle must be fabricated into the forearm when a P.A. is ordered
- Any Hosmer 1-3 / 4" or 2" diameter wrist is suitable
- Select the saddle from Hosmer's conventional or electric elbows

NOTE: When ordering a P.A. with Outside Locking Hinges, the laminating straps are not included.


Please supply the following information when ordering the PA:

- Forearm length from elbow axis to end of wrist. Minimum forearm length: 8" (20.32cm)
- Pigment Color (Standard, Kingsley)
- Right or Left Side
- Wrist Style: Use 1-3/4" or 2" Diameter Wrists
- Elbow Style
- Myoelectric or Switch Control
- Power Supply (Modular Battery Pack)
- Charger

Prehension Actuator P/N	Elbow Size	Wrist Diameter
59165	Medium: E-200	2" (5.08cm)
59166	Medium: E-200	1-3/4" (4.45cm)
59167	Medium: NY Electric	2" (5.08cm)
59168	Medium: NY Electric	1-3/4" (4.45cm)
59163	Large: NY Electric or E-400	2" (5.08cm)
59164	Large: NY Electric or E-400	1-3/4" (4.45cm)
59319	E-500 or E-5500 Elbow Disarticulation Forearm	Custom Order

NY Prehension Actuator System Components

Myoelectric Control Package



- For NY Prehension Actuator

57791 Myoelectric Control Kit

Includes:

- 58671 Myoelectric Processor
- 57804 EMG Electrode Assembly (1 Electrode Included)
- 57809 Ground Electrode Assembly (1 Electrode Included)
- 57815 Electrode Button, Medium (3 Buttons Included)
- 57454 Jumper Connector
- 54487 Gain Adjustment Tool

Switch Controls



- For NY Prehension Actuator

57405 Single Position Pull Switch, NY, For Conventional
1/2" Harness Attachment

53036 Single Position/Single Button Push Switch, NY

71940 Single Position/Single Button Push Switch, OCCC

71977 Three Function Pull Switch, NY, (Not Shown)

NY Prehension Actuator (P.A.) Specifications:

- **Myoelectric Control:** One EMG site and one ground required
- **Switch Control:** Single function control, push or pull switch
- **Power Supply:** AA NiMH rechargeable batteries
- **Weight:** 1 lb (445g) w/9-1/2" forearm and 2" CAPP Delrin Wrist
- **Length:** Minimum length for forearm: 8" (20.32cm)
- **Cable Travel:** Maximum cable travel 1-5/8"
- **Pull:** 12 lbs. or 3 rubber bands
MAXIMUM

Modular Battery Packs

Batteries included


- 56891 5V Remote Battery Pack
- 56356 5V Remote Battery Pack, Dual-In-Line
- 54729 5V Battery Pack in Pouch
- 56776 6V Endoskeletal Battery Pack
- 58058 9V Battery, (Requires Wiring Harness 56876)
- 61859 7.2V, Li-Ion Battery Pack

Battery Chargers

- 54726 5V Battery Charger, 115V (Standard)
- 56993 5V or 6V Battery Charger, 115V
- 71617 5V or 6V Battery Charger, 220V (For International Use)
- 57823 9V Battery Charger
- 60949 Timer For Battery Charger
- 61858 7.2V, Li-Ion Battery Charger Set

Michigan Electric Hook

Michigan Electric Hook



- Child size electric hook
- Assists younger patients with limited excursion in operating a terminal device more effectively
- Utilizes a modified 10P or 10X hook

Features:

- Switch or Myoelectric control
- Available as self contained BE prosthesis with Myoelectric control
- May be used in conjunction with the NY Electric Elbow
- Utilizes nitrile rubber or plastisol covered hooks

Michigan Electric Hook Specifications:

- **Myoelectric Control:** One EMG site and one ground required
- **Switch Control:** Single function control, push or pull switch
- **Power Supply:** AA NiMH rechargeable batteries
- **Weight:** 4oz (113g) 10X Electric Hook (nitrile rubber lined fingers): 5oz (142g) 10P Electric Hook (plastisol covered fingers)
- **Wrist:** Any child size wrist with 1/2-20 threads can be used

Michigan Electric Hook, Nitrile rubber Lined Fingers

70790 Michigan Electric Hook, 10X Right

70792 Michigan Electric Hook, 10X Left

Michigan Electric Hook, Plastisol Covered Fingers

70791 Michigan Electric Hook, 10P Right

70793 Michigan Electric Hook, 10P Left

Myoelectric Control Package

57791 Myoelectric Control Kit

Includes: Processor, 1 EMG Electrode, Ground Electrode, Medium Buttons, and Jumper Connector


Switch Controls

Can be used when the Michigan Electric Hook is used alone or in combination with other electric or myoelectric components

NOTE: See System Components for P.A. on page H114

NU-VA Synergetic Prehensor

NU-VA Synergetic Prehensor Terminal Device



- Lightweight, myoelectrically controlled terminal device
- Opens and closes with true proportional speed
- Pinch force up to 25 lbs (11.3 kg)

Features:

- True Proportional Myoelectric control
- Pinch force maximum 25 lbs (11.3kg)
- Quick and easy interchanging of terminal devices with specially designed wrists
- Convenient 9V battery design, operates a full day on an overnight charge
- Low current drain
- Automatic safety breakaway
- Lyre-shaped, nitrile rubber lined fingers



Gray Terminal Device w/ Black palmar and side pad

- 57840 Synergetic Prehensor, Right, Gray
- 57841 Synergetic Prehensor, Left, Gray
- 60983 Synergetic Prehensor, Utah Euro Wrist Conversion, Right, Gray
- 60984 Synergetic Prehensor, Utah Euro Wrist Conversion, Left, Gray
- 60981 Synergetic Prehensor, Euro Wrist Hosmer Control, Right, Gray
- 60982 Synergetic Prehensor, Euro Wrist Hosmer Control, Left, Gray

Beige Terminal Device w/ Black palmar and side pad

- 57840-S Synergetic Prehensor, Right, Beige
- 57841-S Synergetic Prehensor, Left, Beige
- 60983-S Synergetic Prehensor, Utah Euro Wrist Conversion, Right, Beige
- 60984-S Synergetic Prehensor, Utah Euro Wrist Conversion, Left, Beige
- 60981-S Synergetic Prehensor, Euro Wrist Hosmer Control, Right, Beige
- 60982-S Synergetic Prehensor, Euro Wrist Hosmer Control, Left, Beige

NU-VA System Components

Myoelectric Control Processor Assembly



59860 Myoelectric Processor Assembly Package

Includes:

- 59850 Myo Processor, On/Off Switch and Charging Jack
- 57804 EMG Electrode Assembly, (2 Electrodes Included)
- 57809 Ground Electrode Assembly (1 Electrode Included)
- 57815 Electrode Buttons, Medium,(5 Buttons Included)
- 58058 9 Volt Rechargeable Ni-Cad Battery
- 57823 9 Volt Battery Charger
- 54487 Gain Adjustment Tool

Wrists with Electrical Connectors

- For NU-VA Synergetic Prehensor

- 58052 WE-500N Style
- 58708 WE-500N Style w/Laminating Ring
- 58713 FD 2", Friction Disconnect
- 58704 FM-100 Style w/Laminating Ring

Fabrication

- Call Our Central Fabrication Department For Information

58250 BE Prosthesis Fabrication, Using Synergetic Prehensor.

Options

- 60152 Processor and Electrodes Only, (No Charger)
- 57818 Laminating Form, For EMG Electrode
- 57820 Laminating Form, For Ground Electrode
- 59854 Laminating Form, For Prehensor Processor

Electrode Buttons

- 57814 High Profile, .28" High
- 57815 Medium Profile, .19" High
- 57816 Low Profile, .09" High

NU-VA Synergetic Prehensor Specifications


NU-VA Synergetic Prehensor Specifications

NU-VA Synergetic Prehensor Specifications:

- Proportional Myoelectric Control: Two EMG sites, one to open, one to close and a ground electrode
- Weight: 13.23oz (375g) - Synergetic Prehensor only
- Length: 7-3/32" (180mm), from surface of wrist to top of hook tynes
- Power Supply: 9V rechargeable ni-cad battery. Use only compatible 9V battery charger, 57823, supplied w/ Myoelectric Processor Assembly
- Quiescent Current Drain: 25 micro amps
- Safety Breakaway: 30 lbs. at finger tips
- Pinch Force: 0-25 lbs based on 9V constant power supply. May vary due to patient EMG or condition of battery
- Cycles: 1,300 cycles (more than one full day's normal use) tested on one full charge of ni-cad 9V rechargeable battery. One cycle is defined as full open to full close w/5 lbs pinch force and back to full open. These are relative figures which may vary.

Controls for Electric Components

Prehension Actuator or Michigan Hook Myoelectric Control



57791 Myoelectric Control Kit for One EMG Site,
Single Site/Single Function EMG and Ground Electrodes.
Interfaces w/ Prehension Actuator or Michigan Electric

Hook

NY Electric Elbow Switch Controls



71987 Two Function/Two Button Push Switch, OCCC

Two Function, Extension and Flexion, Internally Mounted

56277 Two Function/Two Button Push Switch, NY Two

Functions, Extension and Flexion, Externally Mounted

57412 Two Function Pull Switch, NY

Operates Extension and Flexion From Conventional

1/2" Harness Attachment

71977 3 Function Pull Switch, NY Activate Three Functions

From One Control Switch: Elbow Extension and Flexion

and Operation of P.A. or Michigan Hook.

Prehension Actuator or Michigan Hook Switch Controls



71940 Single Button Push Switch, OCCC

Single Function, Internally Mounted

53036 Single Button Push Switch, NY

Single Function, Externally Mounted

57405 Single Function Pull Switch, NY Conventional

1/2" Harness Attachment

Controls for Electric Components

- Simple plug-in modules
- Either myoelectric or switch control
- Modular design enables prosthetist to customize an electric system and interchange controls for the best possible application

Power Supplies and Chargers

5V Remote Battery Pack



- Small Size Is Ideal For Mounting This Pack Distal To The Socket Or Directly Above The Elbow Turntable.

56891 5V Remote Battery Pack (4 AA NiMH Batteries)

5V Battery Pack w/Batteries



- Most Commonly Used For Cosmetic Build-Ups When Space Is Not Available Within The Socket Walls.

54729 5V Battery Pack w/Batteries, Case, Pouch (4 AA NiMH Batteries)

Dual-in-Line 5V Remote Battery Pack



- The Long, Narrow Configuration Of This Battery Pack Allows It To Be Placed Within The Socket Walls When Space Above The Elbow Turntable Is Not Sufficient.

56356 5V Remote Battery Pack, Dual-in-Line (4 AA NiMH Batteries)

6V Modular Battery Pack



- Because Of Its Circular Design, This Battery Pack Is Easily Placed And Supported On The Modular Pylon.

56776 6V Modular Battery Pack, Endoskeletal (5 AA NiMH Batteries)

Power Supplies and Chargers

Li-Ion Battery System



- Small size is ideal for mounting this pack distal to the socket or directly above the elbow turntable
- Includes charger AC + DC adapter and battery

61583 7.2 Li-Ion Battery System

Replacement Batteries

- 54737 Battery, 1.25V Rechargeable NiMH,
For Use In 5V or 6V Battery Packs
- 58058 Battery, 9V Rechargeable Ni-Cad,
For Use w/Synergetic Prehensor
- 61584 Battery, 7.2V Li-Ion, For Use w/NY Electric Elbow

Power Supplies/Chargers for Switch or Myoelectric Systems

- Commonly available rechargeable Nickel Metal Hydride batteries will power a device for a full day under normal operating conditions
- Nickel Metal Hydride (NiMH) batteries are environmentally safe

Note: Alkaline batteries can be substituted temporarily when normal recharging is not possible

Battery Chargers

- 54726 5V Battery Charger, 115V,
Designed for Hosmer 5V Modular Battery Pack
- 56993 6V Battery Charger, 115V,
Designed for Hosmer 6V Modular Battery Pack
- 71617 5 or 6V Battery Charger, 220V, For International Use,
or Where 220 Voltage is Customary
- 57823 9V Battery Charger, 115V, Specifically Designed for the
NU-VA Synergetic Prehensor
- 61585 7.2 Li-Ion Battery Charger, Specifically Designed for the
Li-Ion Power System
- 60949 Timer for Battery Charger

Myoelectric and Electric Accessories


Myoelectric Accessories

- 57804 EMG Electrode Assembly
- 57809 Ground Electrode Assembly
- 57814 Button, Electrode, High Profile
- 57815 Button, Electrode, Medium Profile
- 57816 Button, Electrode, Low Profile
- 57818 Laminating Form, EMG Electrode
- 57820 Laminating Form, Ground Electrode

Electric Accessories

- 70409 Forearm Jumper Cable
Specialized Cable Required When Using the
P.A. or Michigan Electric Hook in Conjunction w/Electric Elbow
- 58247 On/Off Push Button Switch Cable Assembly
- 71625 European Plug Adaptor Kit
- 57454 Jumper Connector, for Use w/Myo P.A.
- 54735 3 Pin Connector, Female
- 54736 3 Pin Connector, Male
- 71975 5 Pin Connector, Female
- 56224 5 Pin Connector, Male
- 58196 Shipping Plug, Prevents Activation During Shipping
- 54737 Battery, 1.25v Rechargeable Ni-Cad

MyoLab II Analyzer



- Superior training and fitting aid for myoelectric prostheses
- Readings register immediately without using electrode gel
- Biofeedback reduces learning curve for the patient
- Verify electrode placement inside the socket while patient is wearing the prosthesis
- Ideal as a diagnostic aid and for muscle strengthening or relaxation training

MyoLab II Features:

- Two independent channels
- Locate EMG site on muscle groups for myoelectrode placement
- Allows patient to develop separation of signals through visual feedback
- Compact 7.5" x 4" x 1.2" (19cm x 10cm x 3cm)
- Weights 2 lbs., (445g)

- 61040 Myolab II w/Adaptor
- 61039 Myolab II
- 61038 Adaptor for Hosmer Myo System

